

Anti-DLL3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nyj7

Item Number: CM0022649

Introduction

Rabbit polyclonal antibody targeting human DLL3, suitable for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for Notch signaling and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	DLL3 Rabbit pAb
Remarks/Alias	SCDO1; DLL3
Species	Human
Gene ID (Human)	10683
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 520-618 of human DLL3 (NP_058637.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NYJ7
Protein Weight	65 kDa
Shipping	Ice bag