

Anti-DLL4 Rabbit Monoclonal Antibody For Wb, Ip, Fc, Elisa - Q9Nr61

Item Number: CM0004726

Introduction

Anti-DLL4 rabbit monoclonal antibody validated for WB, IP, FC, ELISA. Targets human DLL4, a Notch ligand involved in angiogenesis and neurogenesis. Suitable for vascular biology and Notch signaling research.

[Learn More](#)

Parameter	Value
Product Name	DLL4 Rabbit mAb
Remarks/Alias	AOS6; delta4; hdelta2
Species	Human
Gene ID (Human)	54567
Gene ID	54567
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 27-524 of human DLL4 (NP_061947.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, FC, ELISA
Cross Reactivity	Human
SWISS	Q9NR61
Protein Weight	75kDa
Shipping	Ice bag