

Anti-Dmtf1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y222

Item Number: CM0023121



Introduction

Rabbit polyclonal antibody targeting human DMTF1, a cyclin-D-binding Myb-like transcription factor involved in cell cycle arrest and tumor suppression. Suitable for Western blot and ELISA applications. Ideal for cancer research.

[Learn More](#)

Parameter	Value
Product Name	DMTF1 Rabbit pAb
Remarks/Alias	DMP1, DMTF, MRUL, hDMP1, DMTF1
Species	Human
GeneID	9988
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 561-760 of human DMTF1 (NP_066968.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y222
Protein Weight	84kDa
Shipping	Ice bag