


Anti-Dnaja2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O60884

Item Number: CM0030888



Introduction

Rabbit polyclonal antibody targeting human DNAJA2 (DnaJ homolog subfamily A member 2). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Useful for chaperone and protein folding studies.

[Learn More](#)

Parameter	Value
Product Name	DNAJA2 Rabbit pAb
Remarks/Alias	DJ3, CPR3, DJA2, DNAJ, DNJ3, RDJ2, HIRIP4, PRO3015, DNAJA2
Species	Human
Gene ID (Human)	10294
Gene ID	10294
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 233-412 of human DNAJA2 (NP_005871.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60884
Protein Weight	46kDa
Shipping	Ice bag