

Anti-Dnal4 Rabbit Polyclonal Antibody For Wb, Elisa - 096015

Item Number: CM0018099

Introduction

Rabbit polyclonal antibody targeting human DNAL4, a dynein axonemal light chain involved in ciliary force generation. Validated for Western blot and ELISA; cross-reacts with mouse and rat. Suitable for respiratory cilia and motor protein studies.

[Learn More](#)

Parameter	Value
Product Name	DNAL4 Rabbit pAb
Remarks/Alias	MRMV3; PIG27; DNAL4
Species	Human
GeneID (Human)	10126
GeneID	10126
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-105 of human DNAL4 (NP_005731.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	096015
Protein Weight	12kDa
Shipping	Ice bag