

Anti-Dnase1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P24855

Item Number: CM0014564



Introduction

Rabbit polyclonal antibody against DNASE1 (Deoxyribonuclease I) for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat; ideal for apoptosis, NETosis, and autoimmune disease research.

[Learn More](#)

Parameter	Value
Product Name	DNASE1 Rabbit pAb
Remarks/Alias	DNL1; DRNI; DNASE1
Species	Human
GeneID	1773
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-235 of human DNASE1 (NP_005214.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P24855
Protein Weight	31kDa
Shipping	Ice bag