

Anti-Dnmt3L Polyclonal Antibody For Wb, Elisa - Q9Ujw3

Item Number: CM0013975



Introduction

Rabbit polyclonal antibody targeting human DNMT3L, a key regulator of DNA methylation. Validated for western blotting and ELISA. Cross-reacts with human and mouse. Suitable for epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	DNMT3L Rabbit pAb
Remarks/Alias	DNMT3L
Species	Human
GeneID (Human)	29947
GeneID	29947
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 151-386 of human DNMT3L (NP_787063.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UJW3
Protein Weight	44kDa
Shipping	Ice bag