

Anti-Dnph1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - O43598

Item Number: CM0004103



Introduction

DNPH1 rabbit monoclonal antibody validated for WB, IF/ICC, and ELISA. Detects human and mouse DNPH1, a 19 kDa nucleotide salvage enzyme that hydrolyzes hmdUMP, preventing genomic damage.

[Learn More](#)

Parameter	Value
Product Name	DNPH1 Rabbit mAb
Remarks/Alias	RCL; C6orf108; dj330M21.3; DNPH1
Species	Human
GeneID (Human)	10591
GeneID	10591
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 75-174 of human DNPH1 (O43598).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O43598
Protein Weight	19kDa
Shipping	Ice bag