


Anti-Dopamine Transporter/Slc6A3 Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q01959

Item Number: CM0014814



Introduction

Rabbit polyclonal antibody targeting human Dopamine Transporter/SLC6A3. Validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for neuroscience and transporter studies.

[Learn More](#)

Parameter	Specification
Product Name	Dopamine Transporter/SLC6A3 Rabbit pAb
Remarks/Alias	DAT; DAT1; PKDYS; PKDYS1; Dopamine Transporter/SLC6A3
Species	Human
GeneID (human)	6531
GeneID	6531
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Dopamine Transporter/SLC6A3 (NP_001035.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01959
Protein Weight	68kDa
Shipping	Ice bag