


Anti-Dp1/Tfdp1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q14186

Item Number: CM0008448



Introduction

High-performance anti-DP1/TFDP1 rabbit monoclonal antibody validated for WB, IHC-P, and ELISA. Reacts with human and rat samples. Ideal for cell cycle and transcription research.

[Learn More](#)

Parameter	Value
Product Name	DP1/TFDP1 Rabbit mAb
Remarks/Alias	DP1; DILC; Dp-1; DRTF1; DP1/TFDP1
Species	Human
Gene ID (Human)	7027
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 311-410 of human DP1/TFDP1 (Q14186).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q14186
Protein Weight	45kDa
Shipping	Ice bag