

Anti-Dppa4 Rabbit Polyclonal Antibody For Wb And Elisa - Q7L190

Item Number: CM0017094

Introduction

Rabbit polyclonal antibody against human DPPA4, validated for WB and ELISA. Detects endogenous DPPA4 across human, mouse, and rat species. Suitable for stem cell and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	DPPA4 Rabbit pAb
Remarks/Alias	2410091M23Rik; DPPA4
Species	Human
Gene ID (Human)	55211
Gene ID	55211
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 15-304 of human DPPA4 (NP_060659.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7L190
Protein Weight	34kDa
Shipping	Ice bag