

Anti-Dpyd Rabbit Polyclonal Antibody For Wb, Elisa - Q12882

Item Number: CM0023985

Introduction

Rabbit polyclonal antibody against DPYD, validated in WB and ELISA. Cross-reacts with human and mouse. Suitable for studying pyrimidine degradation pathways and 5-fluorouracil metabolism.

[Learn More](#)

Parameter	Value
Product Name	DPYD Rabbit pAb
Remarks/Alias	DHP; DPD; DYPD; DHPDHASE; DPYD
Species	Human
GeneID (Human)	1806
GeneID	1806
Immunogen	Recombinant fusion protein containing amino acids 1-173 of human DPYD (NP_000101.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q12882
Protein Weight	111kDa
Shipping	Ice bag