

Anti-Drg2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P55039

Item Number: CM0019246

Introduction

Rabbit polyclonal antibody against human DRG2, a GTPase involved in translation. Validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat DRG2. Ideal for studying GTPase activity and RNA binding.

[Learn More](#)

Parameter	Value
Product Name	DRG2 Rabbit pAb
Remarks/Alias	DRG2
Species	Human
Gene ID (Human)	1819
Gene ID	1819
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 150-260 of human DRG2 (NP_001379.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55039
Protein Weight	41kDa
Shipping	Ice bag