


Anti-Drp1 Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, If-P, Elisa - O00429

Item Number: CM0007079



Introduction

Anti-DRP1 rabbit monoclonal antibody validated for WB, IP, IF/ICC, IF-P, ELISA. Cross-reacts with human, mouse, rat, monkey. Targets dynamin-1-like protein, key for mitochondrial fission and neuronal survival. Ideal for cell biology and disease research.

[Learn More](#)

Parameter	Value
Product Name	DRP1 Rabbit mAb
Remarks/Alias	DLP1; DRP1; DVLP; EMPF; OPA5; EMPF1; DYMPLE; HDYNIV
Species	Human
GeneID (Human)	10059
GeneID	10059
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 561-670 of human DRP1 (NP_036192.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	O00429
Protein Weight	82kDa
Shipping	Ice bag