


Anti-Drp1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - O00429

Item Number: CM0018996



Introduction

Rabbit polyclonal antibody targeting human DRP1 (Dynamin-1-like protein) for WB, IF/ICC, IHC-P, ELISA. Recombinant immunogen (aa 437-736). Cross-reacts with mouse, rat. Suitable for mitochondrial fission, apoptosis, and neurology research.

[Learn More](#)

Parameter	Value
Product Name	DRP1 Rabbit pAb
Remarks/Alias	DLP1; DRP1; DVLP; EMPF; OPA5; EMPF1; DYMPLE; HDYNIV
Species	Human
Gene ID (Human)	10059
Gene ID	10059
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 437-736 of human DNM1L (NP_036192.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00429
Protein Weight	82kDa
Shipping	Ice bag