


Anti-Dsc1 Polyclonal Antibody For Wb, Elisa - Q08554

Item Number: CM0018233



Introduction

Rabbit polyclonal antibody against human DSC1 (Desmocollin-1) for Western blot and ELISA. Validated in human, mouse, and rat samples. Ideal for studying desmosome-mediated cell adhesion and epidermal barrier function.

[Learn More](#)

Parameter	Value
Product Name	DSC1 Rabbit pAb
Remarks/Alias	CDHF1; DG2/DG3; DSC1
Species	Human
Gene ID (Human)	1823
Gene ID	1823
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 480-690 of human DSC1 (NP_077739.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08554
Protein Weight	100kDa
Shipping	Ice bag