

Anti-Dusp1/Mkp1 Polyclonal Antibody For Wb, Elisa - P28562

Item Number: CM0027881



Introduction

Rabbit polyclonal antibody targeting human DUSP1/MKP1 (MAP kinase phosphatase-1). Suitable for Western blotting and ELISA; cross-reacts with human, mouse, and rat samples. Ideal for cell signaling and stress response research.

[Learn More](#)

Parameter	Value
Product Name	DUSP1/MKP1 Rabbit pAb
Remarks/Alias	HVH1; MKP1; CL100; MKP-1; PTPN10; DUSP1/MKP1
Species	Human
Gene ID (Human)	1843
Gene ID	1843
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human DUSP1/MKP1 (NP_004408.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28562
Protein Weight	39kDa
Shipping	Ice bag