

Anti-Dync2Li1 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Tcx1

Item Number: CM0027186



Introduction

DYNC2LI1 Rabbit Polyclonal Antibody validated for Western blot (WB) and ELISA; detects human, mouse, rat cytoplasmic dynein 2 light intermediate chain 1; suitable for cilia biology and intraflagellar transport research.

[Learn More](#)

Parameter	Value
Product Name	DYNC2LI1 Rabbit pAb
Remarks/Alias	LIC3; D2LIC; CGI-60; DYNC2LI1
Species	Human
GeneID (Human)	51626
GeneID	51626
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-352 of human DYNC2LI1 (NP_001180393.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TCX1
Protein Weight	40kDa
Shipping	Ice bag