

Anti-Dynll2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Fj2

Item Number: CM0017559



Introduction

Rabbit polyclonal antibody targeting human DYNLL2, validated for WB, IF/ICC, and ELISA. Detects cytoplasmic dynein light chain 2 involved in intracellular transport across Human, Mouse, and Rat samples.

[Learn More](#)

Parameter	Value
Product Name	DYNLL2 Rabbit pAb
Remarks/Alias	Dlc2, DNCL1B, RSPH22, DYNLL2
Species	Human
GeneID (Human)	140735
GeneID	140735
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-89 of human DYNLL2 (NP_542408.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96FJ2
Protein Weight	10kDa
Shipping	Ice bag