


Anti-Dyrk1A Polyclonal Antibody For Wb, Elisa - Q13627

Item Number: CM0012929



Introduction

Rabbit polyclonal antibody against human DYRK1A, validated for WB and ELISA. Recombinant immunogen, UniProt Q13627. Suitable for studying DYRK1A kinase in DNA repair and transcription.

[Learn More](#)

Parameter	Value
Product Name	DYRK1A Rabbit pAb
Remarks/Alias	MNB; DYRK; HP86; MNBH; MRD7; DYRK1; DYRK1A
Species	Human
GeneID (Human)	1859
GeneID	1859
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 624-763 of human DYRK1A (NP_001387.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13627
Protein Weight	86kDa
Shipping	Ice bag