


Anti-Dystonin (Dst) Polyclonal Antibody For Wb, Elisa - Q03001

Item Number: CM0021297



Introduction

Rabbit polyclonal antibody targeting human dystonin (DST) for Western blot and ELISA applications. Confirmed cross-reactivity with mouse and rat; ideal for cytoskeletal, hemidesmosome, and epithelial adhesion research.

[Learn More](#)

Parameter	Value
Product Name	DST Rabbit pAb
Remarks/Alias	DT, BPA, DMH, EBS3, BP240, BPAG1, EBSB2, HSAN6, MACF2, CATX15, CATX-15, D6S1101, DST
Species	Human
GeneID (Human)	667
GeneID	667
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human DST (NP_001714.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q03001
Protein Weight	861kDa
Shipping	Ice bag