


Anti-E-Cadherin Monoclonal Antibody For Ihc-P, Elisa - P09803

Item Number: CM0010079



Introduction

Mouse monoclonal antibody against E-Cadherin (98kDa, SWISS P09803). Validated for IHC-P and ELISA, cross-reacts with human, mouse, and rat. Ideal for cell adhesion and epithelial research.

[Learn More](#)

Parameter	Value
Product Name	E-Cadherin mouse mAb
Remarks/Alias	Um; UVO; Ecad; ARC-1; E-cad; L-CAM; E-Cadherin
Species	Mouse
GeneID (Human)	999
GeneID	12550
Immunogen	Recombinant protein of mouse E-Cadherin
Source	Mouse
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09803
Protein Weight	98kDa
Shipping	Ice bag