

Anti-Ebp1/Pa2G4 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Uq80

Item Number: CM0010839



Introduction

Rabbit monoclonal antibody targeting human EBP1/PA2G4, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for cell proliferation, ERBB3 signaling, and cancer research studies.

[Learn More](#)

Parameter	Value
Product Name	EBP1/PA2G4 Rabbit mAb
Remarks/Alias	EBP1; HG4-1; ITAF45; p38-2G4; EBP1/PA2G4
Species	Human
GeneID (Human)	5036
GeneID	5036
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-394 of human EBP1/PA2G4 (NP_006182.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UQ80
Protein Weight	44kDa
Shipping	Ice bag