

Anti-Edaradd Rabbit Polyclonal Antibody For Wb And Elisa - Q8Wwz3

Item Number: CM0015681

Introduction

EDARADD Rabbit pAb is a rabbit polyclonal antibody suitable for WB and ELISA, detecting human EDARADD protein. This adapter protein is critical for ectodermal organ morphogenesis and NF-κB signaling.

[Learn More](#)

Parameter	Value
Product Name	EDARADD Rabbit pAb
Remarks/Alias	CR, ED3, EDA3, ECTD11A, ECTD11B, EDARADD
Species	Human
Gene ID (Human)	128178
Gene ID	128178
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-215 of human EDARADD (NP_665860.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8WWZ3
Protein Weight	25kDa
Shipping	Ice bag