

Anti-Edc4 Rabbit Polyclonal Antibody For Wb, Elisa - Q6P2E9

Item Number: CM0018380

Introduction

Rabbit polyclonal antibody targeting human EDC4 (enhancer of mRNA-decapping protein 4). Suitable for WB and ELISA, with cross-reactivity to human and mouse. An essential tool for mRNA decay research.

[Learn More](#)

Parameter	Value
Product Name	EDC4 Rabbit pAb
Remarks/Alias	GE1; Ge-1; RCD8; HEDL5; HEDLS; RCD-8; EDC4
Species	Human
GeneID (human)	23644
GeneID	23644
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 800-1100 of human EDC4 (NP_055144.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6P2E9
Protein Weight	152kDa
Shipping	Ice bag