


Anti-Eef1A2 Polyclonal Antibody For Wb And Elisa - Q05639

Item Number: CM0019289



Introduction

Rabbit polyclonal antibody targeting human EEF1A2, validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying translation elongation and EEF1A2-related neurological and muscle disorders.

[Learn More](#)

Parameter	Value
Product Name	EEF1A2 Rabbit pAb
Remarks/Alias	HS1; STN; EF1A; STNL; DEE33; MRD38; EEF1AL; EIEE33; EF-1-alpha-2; EEF1A2
Species	Human
GeneID (Human)	1917
GeneID	1917
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 184-463 of human EEF1A2 (NP_001949.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05639
Protein Weight	50kDa
Shipping	Ice bag