

Anti-Ehf Rabbit Monoclonal Antibody For Wb, Elisa - Q9Nzc4

Item Number: CM0005406

Introduction

Rabbit monoclonal anti-EHF antibody for Western blot and ELISA applications. Recognizes human ETS homologous factor, a nuclear transcription factor crucial in epithelial cell differentiation and tumor suppression. Useful for cancer research.

[Learn More](#)

Parameter	Value
Product Name	EHF Rabbit mAb
Remarks/Alias	ESE3; ESEJ; ESE3B; EHF
Species	Human
GeneID (human)	26298
GeneID	26298
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of human EHF (NP_036285.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NZC4
Protein Weight	35kDa
Shipping	Ice bag