

Anti-Ehhadh Rabbit Polyclonal Antibody For Wb And Elisa - Q08426

Item Number: CM0030848

Introduction

High-quality rabbit polyclonal antibody against EHHADH (peroxisomal bifunctional enzyme) validated for WB and ELISA. Detects human, mouse, and rat EHHADH (~79 kDa), a key enzyme in peroxisomal fatty acid beta-oxidation.

[Learn More](#)

Parameter	Value
Product Name	EHHADH Rabbit pAb
Remarks/Alias	LBP, ECHD, LBFP, MFE1, PBFE, FRTS3, L-PBE, EHHADH
Species	Human
GeneID (human)	1962
GeneID	1962
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 444-723 of human EHHADH (NP_001957.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08426
Protein Weight	79kDa
Shipping	Ice bag