


Anti-Eif2S2 Rabbit Polyclonal Antibody For Wb, Elisa - P20042

Item Number: CM0024377



Introduction

High-quality rabbit polyclonal antibody against human EIF2S2 (eIF2-beta), suitable for Western blot and ELISA. Cross-reacts with mouse and rat. A reliable tool for investigating translation initiation mechanisms.

[Learn More](#)

Parameter	Value
Product Name	EIF2S2 Rabbit pAb
Remarks/Alias	EIF2, EIF2B, PPP1R67, EIF2beta, eIF-2-beta, EIF2S2
Species	Human
GeneID	8894
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-175 of human EIF2S2 (NP_003899.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20042
Protein Weight	38kDa
Shipping	Ice bag