

Anti-Eif4G2/P97 Polyclonal Antibody For Wb, Elisa - P78344

Item Number: CM0014223



Introduction

Anti-EIF4G2/p97 Rabbit Polyclonal Antibody validated for WB and ELISA in human, mouse, rat samples. Targets translation initiation factor EIF4G2 (DAP5/p97), 102kDa. Ideal for IRES-mediated translation studies.

[Learn More](#)

Parameter	Value
Product Name	EIF4G2/p97 Rabbit pAb
Remarks/Alias	P97; AAG1; DAP5; NAT1; EIF4G2/p97
Species	Human
GeneID (Human)	1982
GeneID	1982
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 750-850 of human EIF4G2/p97 (NP_001409.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78344
Protein Weight	102kDa
Shipping	Ice bag