

Anti-Elongin-A2 (Tceb3B) Polyclonal Antibody For Wb, Elisa - Q8lyf1

Item Number: CM0028203

Introduction

Rabbit polyclonal antibody to Elongin-A2 (TCEB3B) for WB & ELISA. Reacts with human, mouse, rat. For transcription elongation and testis research.

[Learn More](#)

Parameter	Value
Product Name	TCEB3B Rabbit pAb
Remarks/Alias	HsT832; TCEB3B; TCEB3L
Species	Human
Gene ID	51224
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 534-753 of human TCEB3B (NP_057511.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8IYF1
Protein Weight	84kDa
Shipping	Ice bag