


Anti-Emerin/Emd Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P50402

Item Number: CM0017146



Introduction

Rabbit polyclonal antibody against human Emerin/EMD (SWISS P50402), applicable for WB, IHC-P, IF/ICC, and ELISA; detects human and mouse Emerin; ideal for nuclear envelope and Emery-Dreifuss muscular dystrophy studies.

[Learn More](#)

Parameter	Value
Product Name	Emerin/EMD Rabbit pAb
Remarks/Alias	STA; EDMD; LEMD5; Emerin/EMD
Species	Human
GeneID (Human)	2010
GeneID	2010
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human Emerin/Emerin/EMD (NP_000108.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P50402
Protein Weight	29kDa
Shipping	Ice bag