


Anti-Emerin/Emd Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P50402

Item Number: CM0016561



Introduction

Rabbit polyclonal antibody targeting human Emerin (EMD). Validated in WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse. Ideal for nuclear envelope and muscular dystrophy research.

[Learn More](#)

Parameter	Value
Product Name	Emerin/EMD Rabbit pAb
Remarks/Alias	STA; EDMD; LEMD5; Emerin/EMD
Species	Human
GeneID	2010
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human Emerin/EMD (NP_000108.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P50402
Protein Weight	29kDa
Shipping	Ice bag