

Anti-Eml2 Rabbit Polyclonal Antibody For Wb And Elisa - O95834

Item Number: CM0016281



Introduction

Anti-EML2 polyclonal antibody (rabbit) for WB and ELISA applications. Recognizes human EML2 and cross-reacts with mouse. Targets a microtubule-associated protein involved in cytoskeletal regulation.

[Learn More](#)

Parameter	Value
Product Name	EML2 Rabbit pAb
Remarks/Alias	ELP70; EMAP2; EMAP-2; EML2
Species	Human
GeneID (Human)	24139
GeneID	24139
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-415 of human EML2 (NP_036287.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O95834
Protein Weight	71kDa
Shipping	Ice bag