

Anti-Entpd2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y5L3

Item Number: CM0030302

Introduction

Rabbit polyclonal antibody against human ENTPD2 (CD39L1/NTPDase-2) for WB, IF/ICC, ELISA. Targets ectonucleotidase involved in extracellular nucleotide metabolism; suitable for research on cancer, immunology, and cardiovascular disease. - Q9Y5L3

[Learn More](#)

Parameter	Value
Product Name	ENTPD2 Rabbit pAb
Remarks/Alias	CD39L1; NTPDase-2; ENTPD2
Species	Human
GeneID (Human)	954
GeneID	954
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9Y5L3
Protein Weight	54kDa
Shipping	Ice bag