

Anti-Entpd5 Rabbit Polyclonal Antibody For If/Icc And Elisa - 075356

Item Number: CM0028316

Introduction

Rabbit polyclonal antibody targeting human ENTPD5 (CD39L4/NTPDase-5), suitable for IF/ICC and ELISA. Immunogen corresponds to amino acids 149-428; recognizes a 48 kDa protein. Ideal for ER UDPase and glycoprotein folding research.

[Learn More](#)

Parameter	Value
Product Name	ENTPD5 Rabbit pAb
Remarks/Alias	PCPH; CD39L4; NTPDase-5; ENTPD5
Species	Human
GeneID (Human)	957
GeneID	957
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 149-428 of human ENTPD5 (NP_001240.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	075356
Protein Weight	48kDa
Shipping	Ice bag