

Anti-Epb42 Polyclonal Antibody For Wb, Elisa - P16452

Item Number: CM0022946



Introduction

Rabbit polyclonal antibody targeting human EPB42 (Protein 4.2). Validated for WB and ELISA, cross-reacts with mouse. Ideal for research on erythrocyte membrane stability and cytoskeleton.

[Learn More](#)

Parameter	Value
Product Name	EPB42 Rabbit pAb
Remarks/Alias	PA; SPH5; EPB42
Species	Human
GeneID (Human)	2038
GeneID	2038
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-240 of human EPB42 (NP_001107606)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16452
Protein Weight	77kDa
Shipping	Ice bag