


Anti-Ephrin A1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P20827

Item Number: CM0011805



Introduction

Rabbit monoclonal antibody against human Ephrin A1 (B61, LERK-1). Validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse. Ideal for angiogenesis and cancer research applications.

[Learn More](#)

Parameter	Value
Product Name	Ephrin A1 Rabbit mAb
Remarks/Alias	B61; EFL1; GMAN; ECKLG; EPLG1; LERK1; LERK-1; TNFAIP4; Ephrin A1
Species	Human
Gene ID (Human)	1942
Gene ID	1942
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Ephrin A1 (P20827).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P20827
Protein Weight	24kDa
Shipping	Ice bag