


Anti-Epsin 1 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Y6I3

Item Number: CM0003697



Introduction

Rabbit monoclonal antibody against human Epsin 1 (EPN1). Validated for WB and ELISA; cross-reacts with mouse. Detects a 60 kDa band. Suitable for studies on clathrin-mediated endocytosis and membrane trafficking.

[Learn More](#)

Parameter	Value
Product Name	Epsin 1 Rabbit mAb
Remarks/Alias	EPN1; Epsin 1
Species	Human
Gene ID	29924
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 400-500 of human Epsin 1 (NP_001123544.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y6I3
Protein Weight	60kDa
Shipping	Ice bag