


Anti-Erap1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Nz08

Item Number: CM0011592



Introduction

High-quality rabbit monoclonal ERAP1 antibody validated for WB, IHC-P, ELISA. Detects endogenous ERAP1 in human, mouse, rat samples. Suitable for MHC class I antigen processing and blood pressure regulation research.

[Learn More](#)

Parameter	Value
Product Name	ERAP1 Rabbit mAb
Remarks/Alias	ALAP; A-LAP; ARTS1; ERAAP; APPILS; ARTS-1; ERAAP1; PILSAP; PILS-AP; ERAP1
Species	Human
Gene ID (Human)	51752
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 842-941 of human ERAP1 (Q9NZ08).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NZ08
Protein Weight	107kDa
Shipping	Ice bag