

Anti-Ese1/Elf3 Polyclonal Antibody For Wb, Ihc-P, Elisa - P78545

Item Number: CM0026618



Introduction

High-quality rabbit polyclonal antibody against human ESE1/ELF3, validated for Western blot, immunohistochemistry, and ELISA. Suitable for epithelial cell differentiation and tumorigenesis research.

[Learn More](#)

Parameter	Value
Product Name	ESE1/ELF3 Rabbit pAb
Remarks/Alias	ERT; ESX; EPR-1; ESE-1; ESE1/ELF3
Species	Human
Gene ID	1999
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human ESE1/ELF3 (NP_004424.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78545
Protein Weight	41kDa
Shipping	Ice bag