

Anti-Espn Polyclonal Antibody For Wb, Ihc-P, Elisa - B1Ak53

Item Number: CM0020414

Introduction

Rabbit polyclonal antibody against human ESPN, validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Ideal for hearing loss, stereocilia, and actin cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	ESPN Rabbit pAb
Remarks/Alias	USH1M; DFNB36; LP2654; ESPN
Species	Human
GeneID (Human)	83715
GeneID	83715
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 585-854 of human ESPN (NP_113663.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	B1AK53
Protein Weight	92kDa
Shipping	Ice bag