

Anti-Esrrg Rabbit Monoclonal Antibody For Wb, Elisa - P62508

Item Number: CM0005783



Introduction

Recombinant rabbit monoclonal antibody against human ESRRG (Estrogen-related receptor gamma), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for nuclear receptor research.

[Learn More](#)

Parameter	Value
Product Name	ESRRG Rabbit mAb
Remarks/Alias	ERR3; ERRg; NR3B3; ERRgamma; ERR-gamma; ESRRG
Species	Human
GeneID (Human)	2104
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ESRRG (P62508).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P62508
Protein Weight	51kDa
Shipping	Ice bag