


Anti-Estrogen Receptor Beta Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Elisa - Q92731

Item Number: CM0009389



Introduction

Rabbit monoclonal antibody targeting human Estrogen Receptor beta (ESR2), suitable for WB, IF/ICC, IF-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studying hormone signaling, nuclear receptor biology, and cancer.

[Learn More](#)

Parameter	Value
Product Name	Estrogen Receptor beta Rabbit mAb
Remarks/Alias	Erb; ESRB; ODG8; ESTRB; NR3A2; ER-BETA; ESR-BETA
Species	Human
Gene ID (Human)	2100
Gene ID	2100
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-153 of human Estrogen Receptor beta (NP_001428.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92731
Protein Weight	59kDa/54kDa/35kDa/48kDa/42kDa
Shipping	Ice bag