

Anti-Etfb Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P38117

Item Number: CM0022550



Introduction

ETFB rabbit polyclonal antibody validated for WB, IF/ICC, ELISA with reactivity in human, mouse, and rat. Ideal for mitochondrial electron transfer, fatty acid oxidation, and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	ETFB Rabbit pAb
Remarks/Alias	MADD; FP585; ETFB
Species	Human
GeneID	2109
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-230 of human ETFB (NP_001976.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P38117
Protein Weight	28kDa
Shipping	Ice bag