

Anti-Ets2 Polyclonal Antibody For Wb, Elisa - P15036

Item Number: CM0020622



Introduction

ETS2 Rabbit pAb is a polyclonal antibody raised in rabbit against an internal peptide of human ETS2. Validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for transcription research.

[Learn More](#)

Parameter	Value
Product Name	ETS2 Rabbit pAb
Remarks/Alias	ETS2IT1; ETS2
Species	Human
Gene ID (Human)	2114
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 300-400 of human ETS2 (NP_005230.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P15036
Protein Weight	53kDa
Shipping	Ice bag