

Anti-Exosc2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q13868

Item Number: CM0014200



Introduction

Rabbit polyclonal antibody against human EXOSC2, validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Targets exosome complex component involved in RNA processing and degradation.

[Learn More](#)

Parameter	Value
Product Name	EXOSC2 Rabbit pAb
Remarks/Alias	p7; RRP4; SHRF; Rrp4p; hRrp4p; EXOSC2
Species	Human
GeneID (Human)	23404
GeneID	23404
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13868
Protein Weight	33kDa
Shipping	Ice bag