

Anti-Exosc7 Monoclonal Antibody For Wb, Elisa - Q15024

Item Number: CM0003207



Introduction

Rabbit monoclonal antibody against human EXOSC7 (Q15024), suitable for Western blot and ELISA. Detects human and mouse EXOSC7, a non-catalytic subunit of the RNA exosome complex. Ideal for RNA processing and degradation research.

[Learn More](#)

Parameter	Value
Product Name	EXOSC7 Rabbit mAb
Remarks/Alias	p8; EAP1; RRP42; Rrp42p; hRrp42p; EXOSC7
Species	Human
GeneID (Human)	23016
GeneID	23016
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 1-100 of human EXOSC7 (Q15024))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15024
Protein Weight	32kDa
Shipping	Ice bag