

Anti-F4/80 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q61549

Item Number: CM0022110

Introduction

Rabbit polyclonal antibody targeting mouse F4/80 (EMR1), a macrophage-specific adhesion GPCR. Validated for Western blot, immunofluorescence, and ELISA. Demonstrates cross-reactivity with human F4/80. Ideal for immune cell interaction studies.

[Learn More](#)

Parameter	Value
Product Name	F4/80 Rabbit pAb
Remarks/Alias	Emr1; Ly71; F4/80; Gpf480; TM7LN3; DD7A5-7; EGF-TM7
Species	Mouse
Gene ID (Human)	2015
Gene ID	13733
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 367-645 of mouse F4/80 (NP_034260.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q61549
Protein Weight	102kDa
Shipping	Ice bag