


Anti-Fabp7 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O15540

Item Number: CM0017836



Introduction

Rabbit polyclonal antibody targeting human FABP7 (BLBP, B-FABP). Validated for WB, IF/ICC, ELISA. Reacts with human and rat. Detects the 15 kDa brain lipid-binding protein essential for CNS development. Suitable for neurobiology and lipid metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	FABP7 Rabbit pAb
Remarks/Alias	MRG; BLBP; FABPB; B-FABP; FABP7
Species	Human
GeneID (Human)	2173
GeneID	2173
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-132 of human FABP7 (NP_001437.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	O15540
Protein Weight	15kDa
Shipping	Ice bag